
Lumencor SOLA light engine[®] Instruction Manual

Emissions Certifications

This equipment has been tested and found to comply with the limits of EMC directive 2004/108/EC. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications.

Safety Certifications

2006/95/EC

CB Scheme

CE Declaration of Conformity

TUV NRTL Listing

TUV Canadian Listing

TUV European License

Table of Contents

- 1 Introduction**
- 2 Precautions and Warnings**
- 3 Installation and Operating Instructions**
- 4 Spectral Output**
- 5 Routine Maintenance and Trouble Shooting**
- 6 Customer Support**
- 7 Product Specifications**
- 8 Connectors**
- 9 Declaration of Conformance**
- 10 Warranty**

1. Introduction

The Lumencor SOLA light engine (SOLA LE) is designed for laboratory use by bioanalytical researchers and/or developers of life science instrumentation. The light engine provides a white light output by combining 6 discrete, bright, light outputs directly to a sample; or in the case of fluorescence microscopy, to the objective.

Model numbering for SOLA LE is defined as follows: SOLA 6-YYY-ZZ. 6 denotes the number of colors, YYY denotes a unique customer code and ZZ denotes option and revision. The first Z denotes simultaneous output “S” and the second Z denotes the revision level.

The light sources within the SOLA LE are controlled by software; either via a serial interface (RS-232 or USB) to a computer running a Lumencor supplied GUI or a third party microscopy software application. The only manual control is a power switch on the rear panel to turn on/off the 24-30VDC power into the unit. A green power indicator on the rear cover is lit (“ON”) when the power supply is connected to the SOLA LE and the power switch is in the on position.

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A few simple practices will ensure trouble-free operation for the life of the light engine.

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Safety Instructions:

Please read and follow all safety instructions provided **BEFORE** using your new SOLA light engine. Failure to comply with the safety instructions may result in fire, electrical shock, or personal injury and may damage or impair protection provided by equipment. Please save all safety instructions.

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Warning: Statements identify conditions or practices that could result in personal injury.

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Caution: Statements identify conditions or practices that could result in damage to your equipment.

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Safety Items {Mesures de sécurité}:

Warning: DO NOT use an unapproved power supply. The Lumencor-supplied external power supply is recommended for use with the SOLA light engine. Alternate 24-30VDC power supplies may be used provided that the current is limited to 7.9A max. Also, it is imperative that it has output over-current protection, as the power input of the SOLA LE is not fused. The equipment is required to be supplied by a DC power source that has been assessed to meet the requirements of a limited current circuit per clause 9.3 of IEC 61010-1. Connect the AC power cord to a receptacle with a protective safety (earth) ground terminal.

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Warning: DO NOT stare into the output of the light engine. The brightness of this light source is higher than most commercial lighting fixtures and is intended to couple directly into a microscope or other bioanalytical instrument.

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Caution: DO NOT open the unit. There are no serviceable parts inside and opening the light engine enclosure will void the manufacturer's warranty.

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Caution: DO NOT set liquids on the light engine. Spilled liquids may damage your light engine.

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Caution: DO NOT drop the light engine. It contains glass optical components that could be damaged or misaligned by the shock produced by a drop onto a hard surface.

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DISCLAIMER: Lumencor shall not be liable for injury to the user or damage to the product resulting from the SOLA light engine being used in a way for which it was not intended and in complete disregard for all posted safety precautions and warnings.

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3. Installation and Operating Instructions

The SOLA light engine ships with the following list of standard components.

1. SOLA light engine configured with one white output.
2. One 28V/7.9A power supply (Lumencor part no. 27-10003).
3. One 6ft AC power cord for the power supply (for North American customers, Lumencor part no. 29-10002, for UK customers, Lumencor part no. 29-10004 and for European customers, Lumencor part no. 29-10005).
4. One USB-to-RS232 cable (Lumencor part no. 29-10011).
5. One 3mm LLG adapter with one 3mm LLG, 1 meter (Lumencor part no. 82-10012).

When setting the SOLA LE up for use, be sure to place the unit on a hard surface and avoid blocking or restricting airflow at the air inlets or exhaust ports on the enclosure. Restricting the airflow will cause the unit to operate at elevated temperatures and will result in decreased product life and/or premature failure.

Position the unit in such an orientation that allows unrestricted access to the DC power connector. In an emergency, you may need to disconnect power to the unit quickly.

The SOLA light engine can either be controlled by a Lumencor supplied GUI, by 3rd party laboratory software or by the Remote Control Accessory (RCA) via the RS-232 port. The GUI provides a quick and easy way to control your new light engine. You will have the ability to turn the source within the unit on/off, adjust the power of the source from minimum to maximum power. Refer to the photo below of the rear panel for the location of the various connectors.

Rear Panel of SOLA light engine

3.1 GUI Installation

Installation instructions follow for GUI control of the SOLA light engine. Connectivity between the computer and your SOLA LE can be accomplished one of two ways; either using a RS-232 straight-thru cable or using the provided USB-to-RS232 adapter cable. Both methods are covered below.

GUI Installation and Operating Instructions (using RS-232 straight-thru cable)

1. Download the zip file for the SOLA GUI from http://www.lumencor.com/software_control.html.
2. Unzip the file and run setup.exe to install the SOLA GUI.
3. Connect the RS-232 cable between the PC and the SOLA LE.
4. Connect the power supply to the SOLA LE.
5. Toggle the power switch on the rear panel to the ON "1" position. The green LED next to the switch should light.
6. Run the GUI by going to the Program Menu and selecting SOLA Controller.
7. In the COM pulldown menu, select the COM port assigned to serial communications. Generally this is COM1. [If you are unsure which one that is then go to Control Panel, then System, then the Hardware tab. Select Device Manager to see the hardware profile. Expand "PORTS (COM & LPT)" to see which COM port is assigned to the "Communications Port" and select that port in the GUI.] The PC should now have control of all color channels in unison in the SOLA light engine. You will have the ability to turn all colors ON or OFF and adjust the intensity of all sources in unison. Next, hit the INIT button. THIS STEP IS REQUIRED TO RE-ESTABLISH THE COMMUNICATION LINK WHENEVER THE SOLA LIGHT ENGINE IS POWER CYCLED OR WHENEVER THE GUI IS STARTED.


GUI Installation and Operating Instructions (using USB-to-RS232 cable)

1. Download the zip file for the SOLA GUI from http://www.lumencor.com/software_control.html.
2. Unzip the file and run setup.exe to install the SOLA GUI.
3. Connect the USB-to-RS-232 cable between the PC and the SOLA LE. If the PC does not recognize the FTDI device in the cable then you must run the FTDI driver installer, CDM20814_Setup.exe.
4. Connect the power supply to the SOLA LE.
5. Toggle the power switch on the rear panel to the ON "1" position. The green LED next to the switch should light.
6. Run the GUI by going to the Program Menu and selecting SOLA Controller.
7. In the COM pulldown menu, select the COM port assigned to USB-to-Serial communications. [If you are unsure which one is correct then go to Control Panel, then System, then the Hardware tab. Select Device Manager to see the hardware profile. Expand "PORTS (COM & LPT)" to see which COM port is assigned to the "USB-to-Serial Comm Port" and select that port in the GUI.] The PC should now have control of all color channels in unison in the SOLA light engine. You will have the ability to turn all colors ON or OFF and adjust the intensity of all sources in unison. Next, hit the INIT button. (THIS STEP IS REQUIRED TO RE-ESTABLISH THE COMMUNICATION LINK WHENEVER THE SOLA LIGHT ENGINE IS POWER CYCLED OR WHENEVER THE GUI IS STARTED.)

4. Spectral Output

The spectra of the SOLA light engine's white output (BLUE) is shown below along with a 75W Xenon trace (GOLD) and a 120W metal halide lamp output (GREEN).

5. Routine Maintenance and Trouble Shooting

Remove any built-up dust or accumulation on the air intake ports. A vacuum may be used to remove debris so that a steady supply of air is available for cooling. It is recommended that these dust-filters be cleaned by a gentle suction device at least every 6 months and more often in dusty or smoke-filled environments.

There are no user-replaceable components or sub-assemblies in SOLA light engines.

6. Customer Support

T: 503.505.6985

E: les.decker@lumencor.com

W: http://www.lumencor.com/support/software_control

| Specification | Detail |
|---------------------------------|---|
| Temperature | |
| Operating | 32 to 95° F (0 to 35° C) |
| Non-operating | -4 to 158° F (-20 to 70° C) |
| Humidity | |
| Operating and non-operating | 0 to 80% relative humidity, non-condensing |
| Altitude | |
| Operating | 0 to 10,000 feet (3,048 meters) |
| Non-operating | 0 to 20,000 feet (6,096 meters) |
| Dimensions (LxWxH) | 11.0 x 7.0 x 4.0 in / 27.9 x 17.8 x 10.2 cm |
| Weight | 9 lb / 4.1 kg |
| Lifetime | > 15,000 hr |
| Input Power Requirements | 24-30 VDC / 7.9A |
| Warm-up Period | Less than 1 minute |
| Protection | IP Rating of X0 |
| Sound Level | Sound Level at 1 meter < 65db(A) |
| Connectivity | RS-232, USB |
| Warranty | 36 months parts and labor for end users |

M: Lumencor, Inc., 14964 NW Greenbrier Parkway, Beaverton, OR 97006 U.S.A.

7. Product Specifications

The SOLA must be operated and stored within the environmental conditions specified.

8. Connectors

8.2 RS-232 Connector

This port conforms to standard RS-232 Interface Protocol.

| Pins | Definition | DC Characteristics |
|---------------------|-----------------------------------|---|
| 1, 2, 3, 4, 6, 7, 8 | DCD, RXD, TXD, DTR, D3R, RTS, CTS | VCC = 5.0V, V _{low} (max) = 0.8V, V _{high} (min) = 2.4V, I _{low} = 0.5mA, I _{high} = 1.0μA |
| 5 | Gnd | |
| 9 | N/C | |

DB9 Connector Pin Definitions

9. Declaration of Conformance

Manufacturer: Lumencor, Inc.
14964 NW Greenbrier Parkway, Beaverton, OR 97006 USA

We declare under our sole responsibility that the SOLA Light Engine conforms to the following directives and norms:

EMC directive 2004/108/EC
2006/95/EC
CB Scheme
CE Declaration of Conformity
TUV NRTL Listing
TUV Canadian Listing
TUV European License

10. Warranty

The SOLA light engine is backed by a 36 month warranty to end users.