

X-Cite illuminators have been qualified for the following microscopes:

| Microscope Model                  | Light Guide Adaptor<br>(for X-Cite 120 series, 200DC, exacte,<br>XLED1) | Direct Coupling Flange<br>(for X-Cite 120LED) |
|-----------------------------------|---|---|
| <b>CARL ZEISS</b>                 |   |   |
| AxioExaminer                      | 810-00022X  | 810-10022X                                    |
| AxioImager                        | 810-00022X  | 810-10022X                                    |
| AxioObserver                      | 810-00022X  | 810-10022X                                    |
| AxioScope A1                      | 810-00022X  | 810-10022X                                    |
| Axiophot                          | 810-00022X  | 810-10022X                                    |
| Axioplan 1, 2                     | 810-00022X  | 810-10022X                                    |
| Axioskop 1, 20, 50                | 810-00022X  | 810-10022X                                    |
| Axioskop 2, 2FS, 40               | 810-00022X / 810-00035X†  | 810-10022X                                    |
| AxioTech                          | 810-00022X  | 810-10022X                                    |
| Axiovert A1                       | 810-00022X  | 810-10022X                                    |
| Axiovert 10, 35, 100, 135, S100TV | 810-00035X  | TBD   |
| Axiovert 200                      | 810-00022X / 810-00035X†  | 810-10022X                                    |
| Axiovert 400                      | 810-00035X  | TBD   |
| Axiovert 25, 40                   | 810-00037   | TBD   |
| AxioZoom                          | No adapter req'd*   | N/A   |
| Discovery.v8 / .v12 / .v20        | No adapter req'd*   | N/A   |
| LSM                               | No adapter req'd*   | N/A   |
| Stemi SV11                        | No adapter req'd*   | N/A   |
| Stereo Lumar                      | 810-00022X  | 810-10022X                                    |
| <b>LEICA</b>                      |   |   |
| DM2000, DM6000, DMI4000           | 810-00033X  | 810-10033X                                    |
| DMIL                              | 810-00033X  | 810-10033X                                    |
| DMIRB, DMIRM, DMIRE2              | 810-00033X  | 810-10033X                                    |
| DMLB                              | 810-00033X  | 810-10033X                                    |
| DMLFS                             | 810-00033X  | 810-10033X                                    |
| DMR, DMRXA, DMRXA2                | 810-00033X  | 810-10033X                                    |
| MZFLIII, MZ16F, MZ16FA            | 810-00038X  | TBD   |
| MZ10F, M165FC, M205FA             | 810-00038X  | TBD   |
| M10, MZ9.5                        | 810-00038X  | TBD   |
| <b>MOTIC</b>                      |   |   |
| AE31                              | 810-00055   | N/A   |
| BA400                             | 810-00055   | N/A   |

| NIKON                        |                            |            |
|------------------------------|----------------------------|------------|
| 50i, 80i, 90i                | 810-00030X                 | 810-10030X |
| AZ100                        | No adapter req'd*          | N/A        |
| Diaphot, Labophot, TMF       | 810-00030X*                | TBD        |
| E400, E600, E800             | 810-00030X                 | 810-10030X |
| FN-1                         | 810-00030X                 | 810-10030X |
| SMZ1500, SMZ800              | 810-00030X                 | 810-10030X |
| TE2000, all other "TE"       | 810-00030X                 | 810-10030X |
| Ti, Ni, Ci                   | 810-00030X                 | 810-10030X |
| TS100                        | 810-00030X*                | TBD        |
| For "white light TIRF"       | 810-00045                  | TBD        |
| most other Nikon microscopes | 810-00030X                 | 810-10030X |
| OLYMPUS                      |                            |            |
| AX70                         | 810-00032X                 | TBD        |
| BX series: BX50, BX60        | 810-00032X                 | 810-10032X |
| BX2 series: BX41, BX51, BX61 | 810-00032X                 | 810-10032X |
| BX3 series: BX43, BX46, BX53 | 810-00032X                 | 810-10032X |
| CKX41                        | 810-00032X                 | TBD        |
| CK40                         | 810-00032X                 | TBD        |
| CX41                         | 810-00077                  | TBD        |
| IX series: IX50, IX70        | 810-00031X                 | 810-10031X |
| IX2 series: IX51, IX71, IX81 | 810-00031X                 | 810-10031X |
| IX3 series (w/ RFA)          | 810-00031X / 810-00046X †† | 810-10031X |
| IX3 series (w/ RFA-L)        | 810-00031X / 810-00046X †† | 810-10032X |
| MVX10                        | 810-00046X                 | 810-10031X |
| SZX7, 10, 12, 16             | 810-00032X                 | TBD        |

Notes:

\* Additional fittings and/or specific configuration from microscope manufacturer may be required.

† With external devices (e.g. shutters, filter wheels, etc) installed between the adaptor and microscope, use 810-00035X.

†† If using low power objectives (4X, 10X) and viewing via eyepieces, use 810-00046X.